

EmPower™ & MetriCide™ OPA

A United Front

EmPower dual enzymatic detergent and MetriCide OPA Plus high-level disinfectant are the ideal combination for cleaning and high level disinfection of medical instruments.



EmPower

EmPower™ is used to pre-clean critical and semi-critical medical devices prior to terminal sterilization or high level disinfection. It contains two protease enzymes that offer the broadest cleaning action on a variety of protein soils. Inclusion of surfactants provides additional cleaning on carbohydrates, lipid and protein soils. Works fast to remove blood, tissue, mucous and other protein-rich bodily fluids from instruments.

- Dual-enzymatic detergent cleaner
- Less corrosive than distilled water
- Fragrance free option for respiratory-sensitive users
- Gentle on instruments
- Low-foaming formula is ideal for use in automated machines

MetriCide OPA

MetriCide OPA Plus is a high level disinfectant for reprocessing heat sensitive semi-critical medical devices when sterilization is not suitable, and when used according to the Directions for Use. The active ingredient is ortho-Phthalaldehyde and can be used manually or in an automated endoscope reprocessor (AER).

- When used manually, it disinfects in 12 minutes at 20°C and can be reused for up to 30 days*
- When used in an AER, it disinfects in 5 minutes at 25°C in an AER and can be reused for up to 14 days*
- MetriCide OPA Plus Solution Test Strips are fast – see results in only 60 seconds

Ordering Information

Part No.	Description	Packaging
10-4100	EmPower - 1 gallon bottle	4 bottles/case
10-6000	MetriCide OPA Plus - 1 gallon bottle	4 bottles/case
10-602	MetriCide OPA Plus Solution Test Strips -100 strips /bottle	2 bottles/case

*When monitored according to label instructions.
©2019 Metrex Research, LLC. MKT-19-0376