









Our CoolJect mist is an FDA cleared device that offers a cost-efficient, fast-acting, one-handed topical anesthetic application, ergonomically paired with syringe injections.

CJSS-001-REV3 020922







COOLJECTI topical Anesthetic Constacation and Constant and Constant Acceptant

2oz CoolJect Topical Anesthetic Canister

3 Pack - VM02000-03 12 Pack - VM02000-12 (Requires CoolJect Nozzle)

To Start Using CoolJect:



CoolJect Syringe Holder 1-3cc High-Efficiency Nozzle

1 Piece - VM01000 Starter Kit - VM01500-01 (3 Pack and Reusable SH Nozzle)



CoolJect High-Efficiency Nozzle Accessory

1 Piece - VM01100 Starter Kit - VM01500-02 (3 Pack and Reusable Nozzle)

Less Waste & More Shelf Space



Competitor Cans


Competitor cans are twice as large for the same number of uses

Different spray styles require stocking additional canister SKU's



CoolJect Cans

Smaller can size for the same number of uses No need for different spray styles and stocking additional SKU's





Elevated Comfort

Providing the best experience for patients by reducing injection pain.



Targeted Application

Targeted blanched area for immediate, accurate needle penetration.



One-Handed Use

The only device to offer single-handed ability to numb and inject quickly.



Cost-Efficient

CoolJect provides over twice the amount of uses per ounce than the competition.



Fast-Acting

CoolJect mist has the quickest numbing onset & time to injection only takes 4-10 seconds.



Non-Toxic

A safer alternative to ethyl chloride, is hazmat free, and non-flammable.

