


WHAT TYPE OF DOPPLER ARE YOU LOOKING FOR?

Fetal

Calling Point for Fetal Dopplers


CHOOSE PORTABILITY OPTIONS



Vascular

Calling Points for Vascular Dopplers

CHOOSE PORTABILITY OPTIONS



- Fetal Probes**
- Great sensitivity
 - Smooth shape, comfortable for patient, easy to hold and clean
 - D3 – 3 MHz - wide angle, ideal for early fetal hearts
 - D3W - a unique 3 MHz waterproof probe
 - D2 – 2MHz - deeper penetration
 - D2W – 2 MHz - waterproof

- Vascular Probes**
- Great sensitivity
 - Ergonomic shape easy to hold, place in sterile glove finger for use in sterile fields
 - D8 – 8 MHz - best for shallow vessels like radial, ankle, foot, digits, penile
 - D5 – 5 MHz - wider beam, deeper penetration. Good for brachial, carotid, femoral, deep veins