

# Take listening to a whole new level.

3M<sup>™</sup> Littmann<sup>®</sup> CORE Digital Stethoscope



#### Hear more. Way more.

Up to 40x amplification\* and active noise cancellation to reduce background noise.

\*40x amplification based on comparison of analog versus electronically assisted amplified listening with maximum volume at the peak frequency (125 Hz).



Snap tight, soft sealing eartips conform to ears for an excellent acoustic seal and comfortable fit.

## Do more than ever before.

#### Capture

Save and annotate 15, 30, 60 or 120-second recordings in the secure dashboard.

#### Analyze

Create a library of heart sounds to monitor the progression of a disease or for education.

#### **Share**

Securely share recordings with colleagues.

### Hear it the way you want to.

Quickly toggle between analog and amplified listening modes—and use the tunable diaphragms on both sides of the double-sided chestpiece to hear high and low frequency sounds.

## Discover new details with digital.

Connects to Eko software, allowing you to visualize heart sound waveforms in real time. Plus, use AI to aid in the detection of heart murmurs. (Subscription required\*\*. Learn more at ekohealth.com/ai)

Upgrade and listen like never before at Littmann.com/CORE



3M Health Care
3M Medical Solutions Division
2510 Conway Ave., St. Paul, MN 55144 USA
1-800-228-3957 | www.littmann.com

Smart device not included.
\*\*Available in US only

©2020 3M. All rights reserved. 3M, Littmann, the L Littmann logo, the L, and Cardiology IV are marks and/or registered marks of 3M. Unauthorized use prohibited. Used under license in Canada.

Eko trademarks, including CORE, are marks and/or registered marks of Eko Devices Inc.