

Vital Signs® labor and delivery care consumables

## Provide excellent care throughout the birthing process

Vital Signs accessories and supplies for GE Corometrics<sup>™</sup> monitors offer innovation and quality that help clinicians provide exceptional care to expectant mothers and their babies, from preterm monitoring to delivery to discharge. Our antepartum and intrapartum fetal monitoring products provide essential information that helps facilitate quick and effective decision making. Providing quick access to reliable data, our products support the high standard of obstetrics care you and your patients expect.

Q

1

1

1

1

#### External monitoring: toco and ultrasound transducers and belts

Each monitor within the GE Corometrics portfolio has been specifically designed and tested to perform best when used with validated transducers. Our toco and ultrasound transducers come in a variety of options. All transducers come with a two-year warranty

Additionally, our belts are designed to ensure a secure connection of the transducer around expectant mothers, reducing the need to frequently readjust.

Cat. no.

5700HAX

5700LAX

2264HAX

2264LAX

Description

Button style, 8'

Loop style, 8'

Button style, 8'

Loop style, 8'

nd belts	 Cat. no.	Description	Qty. (cs.)
d and	Belts for use with transducer, disposable		
asound warranty.	4425FAO	Loop style with Velcro®, pink/blue	50
ucer			
)ty. <i>(cs.)</i>	4425AAO	Loop style, mesh material, white	10
	4425CAO	Loop style with Velcro, white	10
	8024AAO	Loop style, foam with Velcro ends	150
	2015919-001	Button style, blue/pink-striped	100
	2015827-001	Button style, blue	10

# Internal monitoring: fetal spiral electrodes (FSE) and intrauterine pressure catheter cable (IUPC)

Vital Signs fetal spiral electrodes and intrauterine pressure cables are designed and tested for use with GE Corometrics monitors. When internal monitoring is required, fetal spiral electrodes help monitor the fetal heart rate (FHR) through electric transducers that connect directly to the fetal scalp. FHR monitoring provides the crucial basis for evaluating fetus wellness.

#### Our fetal spiral electrodes design includes:

- Wiring clasp at top of electrode, which allows for steady, one-handed placement and release
- Two-and-a-half-year shelf life
- Sterile, single-use and not made with natural rubber latex

Our IUPC cable helps provide accurate readings of uterine activity.

#### These unique features help enhance patient care:

- Water-resistant connector helps mitigate the risk of erroneous readings.
- Exclusive cable check helps ensure a stable connection.

## Maternal ECG cables and leadwires

Maternal ECG cables and leadwires are available for use with GE Corometrics monitors. Cables come in round or rectangular connection options, and leadwire sets provide three-lead monitoring.

	Cat. no.	Description	Qty. (cs.)	
	FSE and accessories			
	7000AAO	Fetal spiral electrodes	50	
BORERO 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2 2	2464AAO	Attachment pads for FSE cable	50	
	1590AAO	Cable, FSE, round connector to the monitor, 8'	1	
	1591AAO	Cable, FSE, rectangle connector to the monitor, 8'	1	
	IUPC cable			
	1336DAO	Cable, connects IUPC to the monitor, 10'	1	

 Cat. no.	Description	Qty. (cs.)
Y adapter cable for maternal/fetal ECG		
1442AAO	Rectangular connector, 1'	1
Multi-Link <sup>™</sup> E0	CG Cable and Leadwire Kit, 3-Lead, AHA	
1554AAO	ECG cable with round connector, 8' Leadwire, 3-lead, Grabbers, 29"	1



ECG cable with rectangular connector, 8' Leadwire, 3-lead, Grabbers, 29"

## Fetal chart paper

Vital Signs paper products are designed and tested for use with GE Corometrics monitors. With a long history of quality and longevity, our maternal/fetal monitoring paper will provide an accurate record of the delivery that will maintain its integrity for years to come.

	Cat. no.	Description	Qty. (cs.)
n n n n n n n n n n n n n n n n n n n	Fetal chart paper		
	2009828-CAO	30–240 beats per minute (BPM), Z-fold designed for use with Coro 12X, 170 and 250	40
Τότο κάτα (μή τρ. γ 10), 2100 Δ.   Τότο κάτα (μή τρ. γ 10), 2100 Δ.	2009828-DAO	50–210 BPM, Z-fold designed for use with Coro 12X, 170 and 250.	40
	2009828-FAO	50–210 BPM, Z-fold, minute marks designed for use with Coro 12X, 170 and 250	40
	2009828-AAO	30–240 BPM roll designed for use with Coro 14X	6
	2009828-28B	Non-stress test (NST) patient recording forms for fetal charts	150
	2009828-28A	Antepartum nuchal translucency (NT) patient recording forms for fetal charts	300
14	Labels		
	2009828-491	Chart guard labels for records	200

## Other accessories

Our fetal/maternal monitoring accessories have been specifically designed and tested to perform with GE Corometrics monitors. Our products will provide the highest level of quality, performance and safety when used together as a system with Corometrics monitors. Use the remote event marker accessory for simple and reliable annotation.



To learn more, contact your Operative Care sales representative, call 833.327.3284.

#### GLOBAL HEADQUARTERS

Vyaire Medical, Inc. 26125 North Riverwoods Blvd Mettawa, IL 60045, USA

#### Vyaire Medical Oy FI-00510 Helsinki, Finland

## vyaire.com

Romer For U.S., Canada, Argentina, Bolivia, Brazil, Chile, Colombia, Costa Rica, Ecuador, Guatemala, Mexico, Nicaragua, Panama, Peru, Venezuela distribution. Trademarks are the property of their respective owners.

© 2018 Vyaire. Vyaire and the Vyaire logo are trademarks or registered trademarks of Vyaire Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EEC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. CF/VIS/16/0013 (1) (1018)

